



- Foot and hand propel in tilted position
- Rugged construction ideal for heavy duty use
- Shock absorbers for maximum comfort and protection
- Adjustable seat to floor height
- Maneuverable short wheelbase
- Transit-ready option



WHEELCHAIRS FOR FUNCTION AND INDEPENDENCE

www.pdgmobility.com

BENTLEY

THE RIGHT MANUAL TILT-IN-SPACE WHEELCHAIR FOR YOUR NEXT CLINICAL APPLICATION.



Blue Shimmer Tracer

Black Shimmar

Burnt Orange Forest Green Misty Purple











TREND-SETTING FRAME COLORS. YOU CAN ORDER CUSTOM COLORS TOO. CONTACT US FOR DETAILS.







Ideal for hand propulsion



Adjustable seat height, seat depth, back angle and rear wheel



Shock absorbers are for maximum comfort



Long frame option for maximum stability. Protective padding option for additional safety.

Seat width, standard 23" – 32"
Seat width, custom 23" – 32"
Seat depth (adjustable) 16" – 22"
Seat height (adjustable) 90° – 120°
Back height 18" – 26"
Dynamic tilt 20°

Overall width

- for 12" & 16" rear wheels, use seat width plus 9.5"

- for 20" & 24" rear wheels, use seat width plus 10.5" Overall length (no front rigging) Weight capacity Chair weight Transit option available

28.5" - 31"

250 lbs, custom to 450 lbs

39 lbs (without front rigging, wheels, anti-tippers, armrests)

PDG's exclusive TGIF warranty (Thank Goodness It Fits) Lifetime warranty on the frame

PDG recommends all wheelchairs be fitted to the consumer under the guidance of a health care professional knowledgeable of the particular needs of the consumer and the limitations of the product.

USA

PDG Product Design Group PMB # 173 - 1124 Fir Avenue Blaine, WA 98230, USA

CANADA

PDG Product Design Group Unit 102 - 366 East Kent Avenue South Vancouver BC V5X 4N6 T: +1 604 323 9220 F: +1 604 323 9097

Toll Free (in N. America): 888-858-4422 Toll Free Fax (in N. America): 877-323-9093 www.pdgmobility.com